

# Aidal

Homecare Suction Unit

*Make it simple!*

AIDAL is the **brand new suction unit** by Air Liquide Medical Systems.

The attention to each single detail makes it different: it is **studied to meet the exact needs of the user** and facilitate the care of the patient.

Every **user's interaction** with AIDAL **is therefore very simple** but even going beyond the different components, creating a full experience which is **user-friendly, effective and comfortable**.

# Aidal


Homecare Suction Unit



[www.airliquidemedicalsystems.com](http://www.airliquidemedicalsystems.com)

## Contatti

AIR LIQUIDE MEDICAL SYSTEMS S.p.A.  
Via dei Prati 62  
25073 Bovezzo (Brescia)  
Tel. +39 030 20159.11  
Fax +39 030 2098329  
[www.airliquidemedicalsystems.com](http://www.airliquidemedicalsystems.com)





# Aidal

Totally **conceived and developed in Italy** by Air Liquide Medical Systems' team of specialists in respiratory medicine, **AIDAL is a very performing suction unit** with high flow and high vacuum pump being the best solution for any kind of users, professional or not. It is produced with **top quality materials** ensuring a **reliable performance and a long lasting usability**.

## CARE INDICATIONS

- Endotracheal and bronchial aspiration
- Tracheostomy suctioning
- Nasopharyngeal and oropharyngeal aspiration

## User friendly interface

- Easy to read vacuum gauge
- On/Off button with light
- Green led indicator lights up when the device is connected to the mains
- 3 leds battery indicator (green, yellow, red) helps the user understand the remaining battery life and an acoustic signal warns that the battery autonomy is almost over and it is necessary to connect the device to the mains



## Effective and comfortable

- Comfortable vacuum knob for setting the vacuum
- The antibacterial filter is directly attached externally to the cap of the jar, ensuring a fast and simple set up of the suction unit.



## TECHNICAL SPECIFICATIONS

Flow	30 lpm
Vacuum	-0.80 bar
Jar	1000 ml
Dimensions	35 x 28 x 14 cm
Weight	3.6 kg
Operation mode	AC and battery
Battery running time	40 min
Battery recharging time	7 h
IP	IP34



## AIDAL HOMECARE CONFIGURATION

- Suction unit
- External power supply
- Tubing
- 2 Antibacterial filters
- Reusable jar
- Fingertip controller
- Carrying bag